



# Chroma Quant Assessment & Extension





Franklin Guevara • Mariya Hryb • Enisha Sehgal • Tata Serebryany



*CHEM E 546 Winter Quarter Project Update*  
*March 18, 2025*




# Background: Plastic Upcycling




**Need a way  
to analyze  
the products**

# Background: Plastic Upcycling




- Gas collected, 10mL CO<sub>2</sub> internal standard
- Solids collected by filter paper, weighed

## Quantification Workflow












# Julia's Creation: Chroma Quant



[github.com/JuliaH/ChromaQuant](https://github.com/JuliaH/ChromaQuant)

# Data Processing


## Raw Data from Instruments

 MBPR053\_02\_GS1\_MS\_SPEC  
 MBPR053\_02\_GS1\_UA\_Comp\_UPP  
 MBPR053\_02\_GS2\_FID\_SPEC  
 MBPR053\_02\_GS2\_TCD\_CSO  
 MBPR053\_02\_GS2\_TCD\_SPEC  
 MBPR053\_02\_LQ1\_FID\_CSO  
 MBPR053\_02\_LQ1\_FID\_SPEC  
 MBPR053\_02\_LQ1\_MS\_SPEC  
 MBPR053\_02\_LQ1\_UA\_Comp\_UPP

## Component Mixture Breakdown

Carbon Number	Aromatics	Linear Alkanes	Branched Alkanes	Cycloalkanes	Alkenes /Alkynes	Other
1	0	0	0	0	0	0
2	0	0	0	0	0	0
3	0	0	0	0	0	0
4	0	0	0	0	0	0
5	0	0.9357	4.9280	1.4413	0.0955	0.0048
6	0	0.2258	1.6310	0.0587	0.8955	0.0002
7	0	0	0.6621	0.0165	0.0967	0
8	0.0015	0	0.1285	0	0.0117	0
9	0.0001	0	0.0105	0.0001	0.0012	0.0038
10	0	0	0.0002	0.0000	0	0
11	0	0	0	0	0	0
12	0	0	0	0	0	0
13	0	0	0	0	0	0
14	0	0	0	0	0	0
15	0	0	0	0	0	0.0007

## Stacked Bar Chart of Components by Carbon Number



Pre-processed with  
Agilent MassHunter  
Quantitative Analysis  
Software

Option: Manual or Chroma Quant

# How consistent is it? Only some anecdotes

## Liquids Analysis

### Maddie's Manual Quantification

ret	name	formula	Area	Mass of liq (mg)
2.4390576	Isobutane	C4H10	0.449115	0.378199775

### Julia's Automatic Quantification

FID RT	FID Area	MS RT	Compound Name	Formula	m_i
2.443	0.4038	2.439057565	Isobutane	C4H10	0.904900608
2.67	2.6863	2.667417153	Butane, 2-methyl-	C5H12	7.244845844

**Build a software to compare the accuracy of ChromaQuant!**  
(and add an ML model extension to ideally improve performance)

6.4305238	Heptane, 3-methyl-	C8H18	3.10103	2.011377394
6.5737635	Cyclohexane, 1,4-dimethyl-	C8H16	1.23536	1.040296748
6.8663	Cyclopentane, 1-ethyl-1-methyl-	C8H16	0.57493	0.484148596


6.341	0.2098	6.332681415	No Match		
6.433	2.5295	6.430323764	Heptane, 3-methyl-	C8H18	17.51241462
6.576	0.8741	6.573763503	Cyclohexane, 1,4-dimethyl-	C8H16	1.624257867
6.809	0.3706	6.807704841	1,5-Heptadien-4-one, 3,3,6-trimethyl-	C10H16O	0.688651145

Identification was good, but integration results seemed inconsistent

# User Stories for Chroma Quant Extension

User Category	Goals	Needs & Wants	Skill Level
CQ developers (Julia)	Compare the output of manual and automatic analyses	A report on the accuracy of the software.	She knows everything about Chroma Quant because she made it
CQ internal users (Rorrer lab)	Use the assessment tool to confirm by-hand calculations	A way to input additional data into assessment tool	Comfortable with GC/MS FID/TCD. Okay with using CQ
CQ external users (other researchers)	Want to see proof that CQ works before they use it in their lab	Easy-to-read figures and graphs assessing the CQ accuracy	Comfortable with GC/MS FID/TCD. No python, UI only.

# Assessment Tool Flowchart






# Closer Look: Components & Their Tests










## 1. CSV Manual Loader Test


- Simulates a manual analysis CSV file
- Tests that it has correct columns
- Raises error if no file is provided




# Closer Look: Components & Their Tests

## 2. ChromaQuant CSV Formatting Test

 MBPR053\_02\_GS1\_MS\_SPEC  
 MBPR053\_02\_GS1\_UA\_Comp\_UPP  
 MBPR053\_02\_GS2\_FID\_SPEC  
 MBPR053\_02\_GS2\_TCD\_CSO  
 MBPR053\_02\_GS2\_TCD\_SPEC  
 MBPR053\_02\_LQ1\_FID\_CSO  
 MBPR053\_02\_LQ1\_FID\_SPEC  
 MBPR053\_02\_LQ1\_MS\_SPEC  
 MBPR053\_02\_LQ1\_UA\_Comp\_UPP



# Closer Look: Components & Their Tests



# Progress so Far...

```
Code Blame 21 lines (19 loc) · 350 Bytes

1  name: chromaquantenv
2  channels:
3    - conda-forge
4  dependencies:
5    - python=3.10
6    - pip
7    - pandas>=2.2.2
8    - molmass>=2024.5.24
9    - numpy>=1.26.4
10   - scipy>=1.13.1
11   - pillow>=10.3.0
12   - plotly>=5.22.0
13   - pubchempy>=1.0.4
14   - rdkit>=2024.3.5
15   - XlsxWriter>=3.2.0
16   - openpyxl>=3.1.5
17
18   - pip:
19     - chemformula~=1.3.1
20     - ttkthemes~=3.2.2
```

**An environment  
for users to create:**

```
conda env create -f
environment.yml -n
chromaquantenv
```

# ChromaQuantExtension is pip Installable

---

```
pip install -e./ChromaQuantExtension
```




```
conda list
```



```
chromaquantextension    1.0.0          pypi_0    pypi
```

# A Preliminary Interface for the User

---



**Upload Product Composition Sheets**

Manually Calculated Product Composition:

Browse

ChromaQuant Calculated Product Composition:


Browse

**User uploads their  
manually  
calculated csv**


**User uploads  
ChromaQuant  
calculated csv**

# Preliminary Results: Liquid Products of 1 rxn

Manual




Chroma Quant



Not insignificant differences  
in each category

# Preliminary Results: Liquid Products of 1 rxn



## Questions to Answer

- How to compare **MANY** data sets?
- What metrics to **assess**?
- What counts as **“success”**?



# Better File Management

# Big Data Analysis

# Implement ML

# Conda Package

[illegible]

# What We Wish We Wknew Wearlier

---

- Resolving merge conflicts

- Git prune




- Unexpected data formatting

- Dependency issues

- Conda forge vs. pip install



**Thank  
you for  
listening!**



# Full Components Chart


Component Name	What It Does	Inputs (Type Info)	Outputs (Type Info)	How It Uses Other Components	Potential Side Effects
Upload Files	Loads manual data files and CQ input files, converting them to Pandas DataFrames	CSV/Excel (strings, ints, floats)	Pandas DataFrames (strings, ints, floats)	-	File format errors (unsupported file types, encoding issues), corrupt/missing data (empty or incomplete files)
Run CQ Analysis	Runs ChromaQuant on manual data to generate a CQ-processed dataset	CSV of manual data (strings, ints, floats)	CSV of CQ output (strings, ints, floats)	Uses Upload Files	CQ execution failures (software crashes, incorrect input format), unexpected output format
Convert Data to Pandas	Converts CQ and manual output files to Pandas Data Frames	CSV file paths	Pandas DataFrames (strings, ints, floats)	Uses Run CQ Analysis	Parsing errors (wrong delimiters, inconsistent column headers), NaN values due to missing data
Preprocess Data	Cleans and normalizes CQ and manual data for easier comparison	Pandas DataFrames (strings, ints, floats)	Cleaned Pandas DataFrames	Uses Convert Data to Pandas	Data loss (incorrect handling of missing values), data misalignment (mismatched indices or column names)
Compare Compound Names	Performs binary classification to check CQ compound name accuracy	Pandas DataFrames (manual & CQ)	Pandas DataFrame with binary classification (T/F or 0/1)	Uses Preprocess Data	False positives/negatives (misclassifications due to inconsistencies in naming conventions)
Compare Reaction Component Mixture	Computes error metrics (MSE, RMSE, $R^2$ ) for reaction component mixture analysis	Pandas DataFrames (manual & CQ)	Pandas DataFrame with error values (floats)	Uses Preprocess Data	Skewed error metrics (outliers heavily affecting MSE), divergent scales (unit mismatch)
Generate Reports & Plots	Summarizes accuracy metrics, creates visualizations	Pandas DataFrames from comparison components	CSV report & plots	Uses Compare Compound Names & Compare Reaction Component Mixture	Plot rendering issues (if input data is missing or poorly structured), incorrect summary statistics due to bad preprocessing
Train ML Model	Trains a Random Forest model if CQ accuracy > 95%	CSV from reports (strings, ints, floats)	ML model for predicting reaction component mixture	Uses Generate Reports & Plots	Overfitting (if dataset is too small or unbalanced), poor generalization (if CQ accuracy is unreliable)

# Background: Analysis Instruments




GC

## Gas Chromatography




MS

## MASS SPECTROMETRY



FID



TCD

